

HOYA Progressive Fixed Corridor Guide

iD MyStyle 2		
When the fitting height is: (in mm)		Software Chooses: (Corridor Length)
From	To	iD MyStyle 2
0.0	14.9	11
15.0	16.9	12
17.0	17.9	13
18.0	18.9	14
19.0	19.9	15
20.0	99.0	16

iD LifeStyle 4		
When the fitting height is: (in mm)		Software Chooses: (Corridor Length)
From	To	iD Lifestyle 3
0.0	14.9	11
15.0	16.9	12
17.0	17.9	13
18.0	18.9	14
19.0	19.9	15
20.0	99.0	16

iD LifeStyle 3		
When the fitting height is: (in mm)		Software Chooses: (Corridor Length)
From	To	iD Lifestyle 3
0.0	15.9	11
16.0	17.9	12
18.0	21.9	13
22.0	99	14

Array 2		
When the fitting height is: (in mm)		Software Chooses: (Corridor Length)
From	To	Array 2
0.0	16.9	11
17.0	20.9	13
21.0	26.9	15
27.0	99	17